






















# Conventional Skin Closure

## Konventioneller Hautverschluss

### NEEDLE-THREAD COMBINATIONS NADEL-FADEN-KOMBINATIONEN

▼ 3/8 circle, reverse cutting needle  
3/8 Kreis, schneidende Nadel

TYPE OF NEEDLE NADEL TYP	MATERIAL MATERIAL	COLOR FARBE	METRIC USP	0.7	1	1.5	2	3	
			LENGTH LÄNGE	6-0	5-0	4-0	3-0	2-0	
DS 12 	Medilene®		45 cm	51 B 07	51 B 1	51 B 15			
	Nylon		45 cm		51 R 1				
			45 cm		51 C 1				
DS 16 	Medilene®		45 cm	53 B 07	53 B 1	53 B 15			
	Nylon		45 cm	53 R 07	53 R 1	53 R 15			
			45 cm	53 C 07	53 C 1	53 C 15			
	Supramid		45 cm			53 T 1	53 T 15		
			75 cm				53 T 151		
	Polyester		75 cm	53 J 071		53 J 151			
			75 cm		53 K 11	53 K 151			
Silk/Seide		75 cm			53 D 11	53 D 151			
DS 18 	Medilene®		45 cm		54 B 1	54 B 15	54 B 2		
			75 cm		54 B 11	54 B 151	54 B 21		
	Nylon		45 cm		54 R 1	54 R 15	54 R 2		
				45 cm			54 S 15	54 S 2	
				75 cm			54 C 151	54 C 21	
	Supramid		45 cm		54 T 1	54 T 15	54 T 2		
			75 cm				54 T 21		
	Polyester		45 cm		54 K 1	54 K 15	54 K 2		
75 cm					54 K 151	54 K 21			
Silk/Seide		75 cm				54 D 21			
Multipack ▼▼▼▼▼▼▼▼	Polyester		45 cm				654 J 2		

STERILE EO  0482





Packaging Unit: 24 Foils/Box  
VPE: 24 Folien/Box

# Conventional Skin Closure


## Konventioneller Hautverschluss

### NEEDLE-THREAD COMBINATIONS NADEL-FADEN-KOMBINATIONEN

▼ 3/8 circle, reverse cutting needle  
3/8 Kreis, schneidende Nadel

TYPE OF NEEDLE NADELTYP	MATERIAL MATERIAL	COLOR FARBE	METRIC USP	1.5	2	3	3.5	4	
			LENGTH LÄNGE	4-0	3-0	2-0	0	1	
DS 24 	Medilene®	■	45 cm		55 B 2	55 B 3			
			75 cm		55 B 21	55 B 31			
	Nylon	■	45 cm			55 R 2	55 R 3		
			75 cm			55 R 21	55 R 31		
	■	45 cm				55 C 2	55 C 3		
		75 cm				55 C 21	55 C 31		
Supramid	■	75 cm			55 T 21	55 T 31			
Polyester	■	75 cm			55 K 21	55 K 31			
Silk/Seide	■	75 cm			55 D 21	55 D 31			
DS 30 	Medilene®	■	75 cm		57 B 21	57 B 31			
	Nylon	■	75 cm			57 C 31	57 C 351		
	Supramid	■	75 cm			57 T 31	57 T 351		
			45 cm				57 U 35		
	Polyester	■	75 cm			57 K 31	57 K 351		
	Silk/Seide	■	75 cm			57 D 31	57 D 351		
DS 40 	Nylon	■	75 cm			58 C 31			
	Polyester	■	75 cm			58 K 31			
Multipack 	Supramid	■	45 cm			458 T 3			

▼ Straight, cutting needle  
Gerade, schneidende Nadel

TYPE OF NEEDLE NADELTYP	MATERIAL MATERIAL	COLOR FARBE	METRIC USP	1	1.5	2	3	3.5
			LENGTH LÄNGE	5-0	4-0	3-0	2-0	0
GS 63 	Medilene®	■	75 cm			69 B 21		
	Silk/Seide	■	75 cm			69 D 21	69 D 31	69 D 315
	Polyester	■	75 cm				69 K 31	









Packaging Unit: 24 Foils/Box  
VPE: 24 Folien/Box

STERILE EO  0482

# Conventional Skin Closure

## Konventioneller Hautverschluss

▼ 1/2 circle, reverse cutting needle  
1/2 Kreis, schneidende Nadel










TYPE OF NEEDLE NADELTYP	MATERIAL MATERIAL	COLOR FARBE	METRIC USP	1	1.5	2	3	3.5
			LENGTH LÄNGE	5-0	4-0	3-0	2-0	0
HS 12 	Polyester		45 cm	41 K 1	41 K 15			
	Supramid		75 cm		42 T 151	42 T 21	42 T 31	
HS 18 	Polyester		75 cm		42 K 151	42 K 21		
	Supramid		75 cm		43 T 151	43 T 21	43 T 31	
HS 23 	Polyester		75 cm			43 K 21	43 K 31	

Packaging Unit: 24 Foils/Box  
VPE: 24 Folien/Box

# Fascia and Subcutaneous Suture

## Fazien- und Subkutan-Naht

● 1/2 circle, round-bodied needle  
1/2 Kreis, Rundkörper-Nadel

TYPE OF NEEDLE NADELTYP	MATERIAL MATERIAL	COLOR FARBE	METRIC USP	1	1.5	2	3	3.5
			LENGTH LÄNGE	5-0	4-0	3-0	2-0	0
HR 17 	Supramid		45 cm			12 T 2		
	Medifil		75 cm			12 N 21	12 N 31	
	Polyester		75 cm		12 K 151	12 K 21		
HR 22 	Supramid		75 cm		13 T 151	13 T 21		
	Polyester		45 cm		13 K 15			
HR 26 	Polyester		75 cm			14 K 21	14 K 31	14 K 351
	Silk/Seide		75 cm					14 D 351

Packaging Unit: 24 Foils/Box  
VPE: 24 Folien/Box

STERILE EO  0482