


Micro-/ Vessel Surgery

Mikro-/ Gefäß-Chirurgie

NEEDLE-THREAD COMBINATIONS NADEL-FADEN-KOMBINATIONEN

Needle dimensions / Nadel-Abmessungen



	TYPE TYP	CIRCLE KREIS	CURVE KURVE	LENGTH LÄNGE	CHORD SEHNE	NEEDLE Ø NADEL-Ø
Round-bodied needle Rundkörper Nadel	DR 5/15	3/8	126°	5.0 mm	4.1 mm	0.15 mm
	DR 6/20	3/8	121°	6.4 mm	5.1 mm	0.20 mm
	DR 8/22	3/8	123°	8.0 mm	6.3 mm	0.22 mm
	DR 10/23	3/8	120°	10.0 mm	7.8 mm	0.23 mm
	DR 13/30	3/8	135°	12.7 mm	9.9 mm	0.30 mm

● 3/8 circle, round-bodied needle (taper point)
3/8 Kreis, Rundkörper-Nadel

TYPE OF NEEDLE NADEL TYP	MATERIAL MATERIAL	COLOR FARBE	METRIC USP	0.2	0.3	0.4	0.5	0.7
			LENGTH LÄNGE	10-0	9-0	8-0	7-0	6-0
DR 5/15 	Nylon	■	20 cm	02 R 02				
DR 6/20 	Nylon	■	20 cm		03 R 03	03 R 04		
2 DR 8/22 	Medilene®	■	75 cm				204 B 051	204 B 071
2 DR 10/23 	Medilene®	■	75 cm				208 B 051	208 B 071
2 DR 13/30 	Medilene®	■	75 cm					205 B 071

STERILE EO 0482

Packaging Unit: 12 Foils/Box
VPE: 12 Folien/Box